


# User – Manual

## Nikon - A1R MP confocal

(11/27/2017)

From: Karl Ferdinand Ziegler



## Contents

1. Turn On/Off Procedure .....	3
Find the Sample .....	3
2.1. Step 1: Use Brightfield (transmitted light without DIC) .....	4
2.2. Step 2: Use transmitted light with DIC (for higher contrast) .....	5
2.3. Step 3: Use Epi Fluorescence .....	6
3. Software (Confocal) .....	7
3.1. Start Software and Pre-setup .....	7
3.2. Setup the acquisition configuration and capture a 2D image .....	8
3.3. Optical Zoom and move the field of interest .....	10
3.4. Acquire a “Z-Stack” Image (3D-lamge) .....	11
3.5. Different views .....	12
3.6. Scale bar .....	13
3.7. Save data .....	13
3.7.1. For Analysis .....	13
3.7.2. To send as email .....	14
4. Multi Photon .....	15
4.1. Single-color .....	15
4.2. Multicolor – Methods .....	16
4.2.1. Simultaneous with single excitation wavelength .....	16
4.2.2. Sequential with multi wavelength .....	16
4.3. Multicolor and Z-Stack .....	18
5. Troubleshooting .....	19

## 1. Turn On/Off Procedure


1. Power supply for Lasers
2. Key switch for laser
3. Power supply control unit
4. Metal halide lamp (Epi Fluorescence)
5. Halogen bulb
6. Power supply motorized stage
7. Power supply for microscope

Turn on Order 1 to 7 is important. For shutting down: Order in exactly opposite direction



## Find the Sample

### 2.1. Step 1: Use Brightfield (transmitted light without DIC)



1. **After** turn on the Microscope (See 1. Turn On/Off Procedure) and don't start the Software
2. **Choose** an objective
3. Make sure objective is clean and **prepared** with Oil/water if necessary
4. **Place** sample on stage
5. Make sure **Filters** are out of light Path [Panel A]
6. **Turn On** the Transmitted light (Button) and change the intensity (Knob) [Panel C]
7. Switch the **Filter turret** to an empty position [Panel D] (dashed line) by using the button [Panel E]
8. Use the **Eye pieces** to find the Sample



### 9. Control Stage (X,Y,Z) see Panel F:

- a. Rotate knob [1] change the moving speed for X,Y movement (current speed shown in Area 1b)
- b. Translate knob [1] move stage in X,Y position respectively
- c. Press button [2] change the moving speed for Z movement (current speed shown in Area 2b)
- d. Rotate knob [3] move stage in Z position

## Comments:

1. **If contrast is not high enough use DIC (see description below)**
2. **To reduce photo bleaching turn the light off if don't needed and use only as much intensity as necessary.**

## 2.2. Step 2: Use transmitted light with DIC (for higher contrast)




1. **Find** sample with bright field (2.1 Step: 1)
2. **Replace** Step 7 in 2.1 with the following 2 Steps
3. Make sure:
  - a. **Polarizer** is in light path [Panel B]
  - b. **Wollaston Prism** is in place [Panel C]
  - c. Under the Objective is another **Wollaston Prism** [Panel D]
  - d. Switch [Panel F] the Filter turret to **Analyzer** [Panel E]
4. Rotate the polarizer to change the DIC **contrast** [Panel G]




### Comments:

1. Some objectives are not available with Wollaston Prism. Both 20x + 60x will have it. If you want DIC for other objectives (e.g. 40x), please ask core facility.

### 2.3. Step 3: Use Epi Fluorescence



1. **Step 1 – 4** from 2.1. (Step 1: Use transmitted light without DIC)
2. Open the **shutter** [Panel B knob I] for the metal halide lamp and **change** intensity [Panel B knob II]
3. Switch Filter turret to one of the Filters (**DAP, FITC or TxRed**) [Panel C, D and E] by pressing the button [Panel D]
4. Use the Eye pieces to find the sample



#### Comments:

1. To reduce photo bleaching turn the light off if don't needed and use only as much intensity as necessary.
2. DIC is not available in this mode.

### 3. Software (Confocal)

After the:


1. Sample is in focus
2. Transmitted light is turned off (reduce photo bleaching).
3. Shutter of the metal halide lamp is closed (reduce photo bleaching).

Now it's ready to start the Software and use confocal.

#### 3.1. Start Software and Pre-setup

##### LUTs:

Look Up Table to see current intensity values and scaling (3.2. Setup the acquisition)



##### General Acquisition settings:

Change scan speed, laser power etc. (3.2. Setup the acquisition)

##### Optical Zoom:


To increase magnification but decrease the field of view (3.3. Optical Zoom)

Two ways to configure the confocal acquisition settings

- a) Select optical configuration [Panel A, Button 1] labeled "A1-488-561-640 TD"

OR


- b) Open a previously acquired image. Right click on the image → Reuse camera settings [Panel B]



### 3.2. Setup the acquisition configuration and capture a 2D image




1. Choose **Light Sources** you want to use by checking/unchecking the boxes next to it [Panel A, Button 1 respectively]:
  - a. Green Bar: Laser 488 nm and Detector
  - b. Red Bar: Laser 561 nm and Detector
  - c. Blue Bar: Laser 638 nm and Detector
  - d. Transmitted light and Detector
2. If **Remove Interlock** [Panel B, Button 1] is red → press it to remove interlock. If the red color does not disappear then make sure L100 on the front of the microscope is selected
3. Press **Scan** [Panel B, Button 2] button to turn on the laser and get a **live image**. Press it again and you turn the laser off and freeze the image.
4. To prevent **pixel saturation** on the detector and achieve a good **signal/noise** ratio the intensity values in the LUTs [Panel C, Field 1 and 2] should be between 2000 and 3000 [Panel C, marked value]. You can achieve that on three different ways:



- a. Adjust the **gain** [Panel A, Area 2 respectively for each la  
Effect: Increase of gain increase signal and noise at the same time
- b. Adjust the **scan speed** [Panel B, Area 8].  
Effect: Reduce speed increase signal but increase photo bleaching (see also the following step 5)
- c. Adjust the **Laser power** [Panel A, Area 4 respectively for each laser]  
Effect: Increase of laser power increases photo bleaching

**(It's necessary to make compromise between these 3 settings above)**

5. Choose the **Images Size** in Pixel [Panel B, Area 9]. This is pixel density and contributes to resolution image quality. It is not possible to use all combinations of scan speed and image size. You have to take that in consideration in Step 4. A standard setting is 512 or 1024.



6. There are three additional settings (Step 4a) that has an effect on the **scan speed**:
  - a. Switch between **Galvano** and **Resonant** scan mode [Panel B, Area 4]  
Effect: Resonant increase speed and noise as well.
  - b. Switch between **uni-** and **bidirectional** scan [Panel B, Area 6]  
Effect: Bidirectional double speed but increase noise
  - c. Change the **Averaging** [Panel B, Area 10] between normal, line and integral. It is only recommended to use Averaging only when capturing image and not during live.  
Effect: Line and integral reduce scan speed but increase signal/noise ratio (clean up but slow)
7. Change the **Pinhole size** (AU) in [Panel B, Area 12]. This value is wavelength dependent. It is recommended to use a value of 1.2 AU (next to the text box)
8. After you are done with all settings and you are ready to acquire an image the press **Capture** [Panel B, Button 13]

## Comments:



1. If you **lost your sample**: the button **Eye Port** [Panel B, Button 3] will turn off the laser and detector and switch to the Eye Port. Now it possible to use the transmitted light or the Epi-Fluorescence. After pressing the button again the old settings get restored.
2. **Unit of scan speed** [Panel B, Area 7]. It is recommend to use Frame/sec.
3. **Shows** the current fps and Frame time [Panel B, Area 11].
4. To reduce **photo bleaching**: freeze the image as often as possible (that will turn of the laser light)
5. There are two ways to adjust image scaling to the monitor [Panel C, Area 1 and 2]
  - a. Keep Auto scale [Panel C, Area 3]
  - b. Auto scale [Panel C, Area 3]

### 3.3. Optical Zoom and move the field of interest

Increase the magnification independent from the objective but decreases field of view as well.





1. **Zoom** in a field of interest [Panel A, area 1] by entering a zoom value or moving the slider [Panel A, area 2] or changing the size of the field with the mouse
2. After you changed the size the green square [Panel A, area 1] becomes red. With a right-click in the area you can **accept** the new field and it turns into green again
3. You can also **move** the field by drag and drop in the entire scan field. But it is recommended to keep the field of interest in the middle of the scan area. Don't move it too far to the corners
4. If you want to **move** the center of your field of interest then make a right-click at this point in the window [Panel C, Area 1] and choose "Move this Point to Center" [Panel B, Button 1]



### 3.4. Acquire a “Z-Stack” Image (3D-Image)


After you found the right acquisition settings (3.2) and the an area of interest (3.3) you can acquire a Z-Stack image (3D Image) by pressing the Button 1 in Panel A. After that the window shown in Panel A, Area 2 or in big, Panel B will appear.




1. Choose a Path to **save** the files and a Filename [Panel B, Area 1]
2. Uncheck all **boxes** [Panel B, Area 2] except for Z.
3. Click on **symbol** left beside Z [Panel B, Area 2] 
4. Now you have three **option** to acquire a Z-Stack image [Panel B, Area 3]: 



- a. **Left** [Panel C]: Move the objective up and down by scrolling mouse wheel or using the stage control unit (see 2.1). When you found the lower limit of your region of interest then press the button “Bottom”. For the upper limit press the button “Top” [Panel C, Area 1]. Choose your step size or your amount of *steps* [Panel C, Area 2]. See comments.



- b. **Middle** [Panel D]: Enter the thickness of your Z-Stack image [Panel D, Field 1]. Your current position is symmetric in the middle of this value. Choose your step size or your amount of *steps* [Panel D, Area 2]. See comments.




- c. **Right** [Panel E]: Enter the thickness of your Z-Stack image by typing a value relative below to your current position and a value relative above of your current position [Panel E, Area 1]. Thus, your current position is asymmetric in the middle of this two values. Choose your step size or your amount of *steps* [Panel E, Area 2]. See comments.

5. Press **Run** [Panel G, Field 3 (next page)] to acquire a Z-Stack image. Please pay attention for the comments on the next page.

## Comments:

1. You can choose the value for Step size or the number of steps. If you enter a step size then the program calculate the amount of steps and opposite. But it is **recommended** to use always as step size the value between the field Step and Steps. You can set this value by clicking the field [Panel F, Field 1]
2. Make always sure that for Z-Device is “Ti ZDrive” chosen [Panel G, Field 1].
3. Make always sure that the **box** beside “Close active Shutter during Z Movement” is checked [Panel G, Area 2].



### 3.5. Different views


<p><b>Slice View:</b> This view displays orthogonal XY, XZ and YZ projections of the image sequence. (Requires Z or T dimension)</p>	<p><b>Volume View:</b> This view creates a 3D model of the acquired object. (Requires Z dimension)</p>	<p><b>Tiled View:</b> This view displays frames of the selected dimension arranged one next to other. (Requires Z, T or XY dimension)</p>

Clicking the buttons open **new** windows with the selected view.

Drag and drop this buttons in the current window **changes** the view.


### 3.6. Scale bar

If you want to **add** a scale bar in your captured image then click on the button on the right side of your acquired image and choose between **Vertical** and **Horizontal** to burn the scale in your image or move the scale bar around.



### 3.7. Save data

#### 3.7.1. For Analysis




1. Two ways to **save** your data:
  - a. **Press** one of the button [Panel A; Area 1]
  - b. **Click** on File [Panel A; Button 2] and then select save as
2. Select always the Image File format **ND2** as shown in Panel B, Area 1.
3. After selecting the Image File Format ND2 press **save** [Panel B; Button 2]

#### Comments:


1. This saves the image in a **raw format** and the microscope configuration as well
2. That allows you to **reuse the camera** settings as described in point 3.1. (Start Software and Pre-Setup) and to change your analysis after acquire an image.
3. The data format is quite too big for sending as **email** attachment. If you want to send images as email then follow the Section 3.7.2. below.



### 3.7.2. To send as email



1. Bring the image you want to send to the **front of the GUI** so it is visible
2. Now you have **two options** to generate an image:
  - a. Option 1:
    - i. Press the Edit Button in the **top left corner** [Panel A; Button 1] to open the drop down menu
    - ii. Select **“Create Full View Snapshot (8bit RGB)”** [Panel A; Area 2] to generate a 8bit RGB image
  - b. Option 2:
    - i. Press the key **combination “Shift + Z”**
3. Now, **select File → Save as** in the top left corner like it is shown in 3.7.1. or press one of the two save symbols. 
4. Select the **Image File Format tiff** [Panel B; Area 1]
5. Make sure that the **boxes are checked** as shown in Panel B; Area 2
6. Press **save** [Panel B; Button 3]




### Comments:

1. The Settings shown in Panel B are only **recommendation**. Feel free to make your own choices.


## 4. Multi Photon

### 4.1. Single-color


**TURN OFF THE ROOM LIGHT AND CLOSE THE DOOR BEFORE YOU USE MULTI PHOTON**



1. To **pre-setup** the system for Multiphoton Microscopy press the Button MP [Panel A; Button 1].
2. After couple of seconds the **System is configured**. An indicator is that the system changed the wavelength from all laser sources to 800 nm [Panel B; Area 1 - 4].
3. Unselect three boxes and leave one [Panel B; Area 5 - 8]. You will have ability to adjust HV and laserpower
4. Changes in the **Laser power** count for each channel. The number "800" on the bottom left side of each channel are a small indicator. [Panel B; Area 9 - 12]. In this mode you can use only one Laser source.



5. After the System is setup you can **open the GUI** (Graphical User Interface) to control the multi photon configuration. Press "right click" somewhere in the black screen and choose "Acquisition Controls" → "A1 MP GUI" [Panel C; Area 1+2].
6. Now the **Window named "A1 MP GUI"** [Panel E] should appear. Sometimes the window appears beside the "A1 Compact GUI" in the top right corner on the screen [Panel D; Area 1]. You can keep it there or take it out per drag and drop.



- To turn the emission **laser on** press the “On” button [Panel E; Area 1]. After that a yellow light right beside the button will indicate the warm up process of the laser [Panel E; Area 1]. After the light turns into green the laser is ready to use.
- To turn the excitation On/Off use the two buttons for the **shutter** [Panel E; Area 2]. The left button (yellow circle) indicates an open shutter and the right button (grey circle) indicates a closed shutter.
- Before you scan the sample, you need to select the optimum wavelength. To **change the excitation** of the laser use the controller name “Wavelength” [Panel E; Area 3]. You can type in values (Unit nm) between 700 - 1000 on the right side or move the slider to adjust the wavelength.
- Make sure Auto is selected for laser alignment [Panel E; Area 4].
- After you setup all other settings, e.g. the scan speed and images size you can scan the sample by pressing the button **scan** or acquire an image by pressing the button **Capture**.

## 4.2. Multicolor – Methods


### 4.2.1. Simultaneous with single excitation wavelength

- Make sure appropriate number of boxes are checked per fluorophore [Panel B]. Allows individual control of HV for each channel.
- Follow procedure for single color PM however allows individual HV control


**Comments:** With this method it is not possible to use the optimum wavelength for each fluorophore. It is necessary to make compromises between the different probes. However this method is the fastest.

### 4.2.2. Sequential with multi wavelength


- Click in the top left corner on “calibration” → “Optical Configuration” or press “Ctrl+N” [Panel A; Area 1] to open the **Optical configurations** [Panel D].
- Click on the left side in the “optical configurations” on MP [Panel B; Area 1] and click on **duplicate** [Panel D; Button 2].
- After that **change the name** of the duplicate [Panel D; Area 3] in “MP-green-(your initials)” for green, “MP-red-(your initials)” for red and so on.
- Then **close** the “Optical configurations”.
- Choose on the top of the screen one of your **generated** “Optical configuration” [Panel B]



- Setup the system for one color according to **point 4.1. Single-color**. (Wavelength, Power, Gain etc.)
- Then **save** by right click on the “Optical configuration” on the top and select “Assign current camera setting” as shown [Panel C; Area 1]. For each fluorophore you are using, repeat the creation of independent optical configuration and all the appropriate laser, detection and scan settings. The red exclamation mark indicate a not saved setting.



C



8. **After you configured** the system for each wavelength click on the left side of the screen on the symbol shown in *Panel E; Area 1*. ND acquisition wizard appears.


9. Within the ND wizard **uncheck all boxes** [*Panel F; Area 1*] except the following symbol:



10. **Add your different "Optical configurations"** on the left side on the window [*Panel G; Area 2*] by checking



11. **Change the order** of the "optical configurations" by clicking the little arrow on the right side in of *Area 2* in *Panel G*.

12. To start the experiment press **run now** [*Panel G; Button 4*]



### 4.3. Multicolor and Z-Stack

1. Setup a Z-Stack experiment according to point 3.4. Acquire a “Z-Stack” Image (3D-Image)
2. Setup the multicolor MP experiment according to point 4.2. Multicolor
3. For faster z-stack acquisition make sure the Symbol for Z-Stack [Panel F] is in the left side of the Multicolor by changing the order [Panel G; Area 3]
4. Press Run Now [Panel G; Button 4]



### Comments:

1. **BEFORE YOU USE THE MULTI PHOTON OPTION ON THIS MICROSCOPE IT IS VERY IMPORTANT TO TURN OFF THE ROOM LIGHT AND CLOSE THE DOOR. OTHERWISE YOU DAMAGE THE DETECTORS.**
2. Before you **leave the room** switch the microscope mode back to confocal in the top of the screen [Panel F; Area 1]. The microscope will change the active detectors.



3. It is **highly recommended** to keep the order for the combination Multicolor and Z-Stack as shown in Point 4.3. otherwise the acquisition takes a very long time because the system have to change the wavelength between each Z-Image.
4. MP recommended objectives are:
  1. Plan Fluor 20x MI
  2. APO Lambda S LWD 40x WI
  3. Plan APO IR 60x WI

## 5. Troubleshooting



If after starting the Software an error message appears and one or more of the device are not recognized by the software then follow the next steps:

1. Turn off the **Software** → Wait 30 sec.
2. Turn off the **control unit** [*Panel A, Field 1*] (Step 3 of the Turn On/Off procedure; Page 3) → Wait 1 min
3. Make sure all **7 power switches** are turned on except the control unit [*Panel A, Field 1*]:
  - a. If that is the case **skip** to step 4
  - b. If that is not the case: turn on **all missing device** except the control unit (Step 3 of the Turn On/Off procedure; Page 3) by following the turn on procedure (see 1. Turn On/Off procedure; Page 3).  
**IMPORTANT: DON'T TURN EVERTHING OFF.**
4. Wait 1 min
5. Turn on the **control unit** [*Panel A, Field 1*] → Wait 1 min
6. Start **Software**

### Comments:

1. **If the problem still appears ask** Andy Schaber (Facility manager) or Ferdinand Ziegler (TA in Facility) for support.
2. If you **switch off the Metal halide lamp** (Step 4 of the Turn on Procedure; Page 3) **PLEASE** don't turn it immediately on again. That can damage the light source. Wait at least 15 min before you turn the lamp on again. You can use the Software without the metal halide lamp but you can't use the Epi fluorescence to find the sample.